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**Micro and Semi-micro-SPALTROHR®-Columns for the gentle and high efficient separation by distillation**

---

**FISCHER**  
technology ... a brand of i-Fischer Engineering® GmbH


Dachdeckerstr. 1 • D-97297 Waldbüttelbrunn  
Tel.: +49-931-46797-0 • Fax: +49-931-46797-29  
e-mail: [fischer-technology@iludest.de](mailto:fischer-technology@iludest.de) • [www.fischer-technology.de](http://www.fischer-technology.de)

### Micro and Semi-Micro-SPALTROHR<sup>®</sup> - Columns for the gentle and highly efficient separation by distillation

Laboratory equipment for the gentle and highly effective separation by distillation

*FISCHER* SPALTROHR<sup>®</sup>-columns rank world wide first amongst the laboratory equipment for separation by distillation. Important features for various Models in an over view:

- Highest separation efficiency and lowest flask pressure compared to all other known types of distillation columns. SPALTROHR<sup>®</sup>-columns permit flask pressures as low as 0.11 mbar.
- Different models permit optimum adjustment to the particular separation problems and charge quantities.
- All-glass construction  
All parts coming into contact with the distillation product are made of glass. Spherical joints are grease-free.
- Minimum hold-up

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**Fig. 1 Axial section of external tube (A-B)**

### **Semi-micro-SPALTROHR<sup>®</sup>-columns HMS 300 C and HMS 500 C**

Semi-micro columns are suitable for charge quantities of approx. 10 to 200 respectively 350 ml.

Corresponding to the larger design, the column is constructed for a higher separation efficiency and automation. A column head with automatic vapour divider serves for controlling the reflux ratio. The requested times of reflux and take-off, in the range of 1:1 to 1000:1 and reciprocal are pre-selected via an electronic time controller. The valve anchor is lifted by means of a solenoid coil.

Additionally to the measurement of the specific column head temperature, an automatic boiling limit can be pre selected. If the actual head temperature exceeds during the distillation the boiling limit, which is set in the control device, the distillate take-off is interrupted automatically. The flask is heated by an oil bath with integrated magnetic stirrer. The liquid in the oil bath and the sample are stirred synchronously via two rotary magnets. This ensures that the evaporation is carried out steadily without boiling retardations

The electric lifting platform enables that the height of the oil bath can be continuously varied according to the requested immersion depth of the flask.

### **FISCHER Distillation Control Device MINITRON [DC/MIC]<sup>2</sup>**

The processor based control and monitoring device combines up-to-date distillation technique with the advantages of modern microelectronics and is distinguished by operational safety, reliability and ease of operation. The safety devices are designed for unattended operation and automatically end the distillation when a failure arises.

**Easy Attendance** by Touch-Screen and Process Data Storage:

Process-/Measuring-Data are appearing on a large 10.1" display (210 x 155 mm) with a resolution of 640 x 480 pixels.

For the handling of different process mixtures up to **100 data sets can be stored.**

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**The control device provides the following functions:**

- start-up mode with increased heat input and under total reflux
- automatic change of pre-distillate cut
- withdrawal of various cuts into separate receivers
- intermitted operation for getting an optimum in distillate quality and -quantity
- automated process termination at adjustable cooling water flow delay

**Operation without Supervision** due to continuous control of:

- cooling water flow
- temperature limits

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
**HMS 500 AC – Fully automatically controlled SPALTROHR<sup>®</sup>- column with fraction collection  
system and data output**

### **HMS 500 AC**

**Fully automatically controlled  
SPALTROHR<sup>®</sup>-column with fraction  
collection system and data output**

Proven supplementation for the preliminary examination of crude-oil samples following ASTM 2892 and for the fine fractionation of boiling cuts.

### **FISCHER Distillation Control Device MINITRON [DC/MIC]<sup>2</sup>**

The processor based control and monitoring device combines up-to-date distillation technique with the advantages of modern microelectronics and is distinguished by operational safety, reliability and ease of operation.

The data presentation software **MINITRON MIC-TOP**, running under WINDOWS<sup>™</sup> on a PC existing on-site: measuring values can be displayed graphically, stored, recalled at any time and printed.

Further a remote control of the distillation unit via RS 232 of the PC is also possible.

The safety devices are designed for unattended operation and automatically end the distillation when a failure arises.

**Easy Attendance** by Touch-Screen and Process Data Storage:

Process-/Measuring-Data are appearing on a large 10.1" display (210 x 155 mm) with a resolution of 640 x 480 pixels.

For the handling of different process mixtures up to **100 data sets can be stored.**

**The control device provides the following functions:**

- start-up mode with increased heat input and under total reflux
- automatic change of pre-distillate cut
- withdrawal of various cuts into separate receivers
- intermitted operation for getting an optimum in distillate quality and -quantity
- automated process termination at adjustable cooling water flow delay

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- cooling water flow
- temperature limits
- level control (intermediate receiver)

For the distillation of volatile substances as well as fraction cuts over a large boiling range or for prolonged distillation, our package includes a distillate fraction take-off system.

It allows to transfer the distillate from vacuum to atmospheric conditions also.

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
**Micro SPALTROHR<sup>®</sup>-column MMS 155 AC**  
**and MMS 255 AC**

The smallest and very compact glass apparatus consists of a flat-bottom flask for processing 0.5...5 ml and a micro-SPALTROHR<sup>®</sup>-column in a silvered high-vacuum jacket. The reflux ratio can be infinitely varied by means of the needle valve. The filling and emptying of the flask is done through the screw cap for the temperature sensor. In case of the larger column MMS 255 AC the flask is separated via a grease-free spherical joint. An oil bath with built-in magnetic stirrer serves as flask heater. Oil bath liquid and flask charge are simultaneously stirred via two rotary magnets. This ensures that evaporation is carried out continuously without boiling retardations.

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### Fraction Collector

An easily operated unit. The revolving insert holding up to 22 collecting tubes can be removed by operation of the quick release spring and kneeling of the lower part of the collector. The collecting tubes can now be exchanged. The top part of the change of tubes is effected by turning the wheel at the base part. A sample can be taken from any of the fractions through an injection part, also under vacuum.

### Sample Changer

This practical device enables the operator to change single collecting tubes also under vacuum. The distillate is collected via a greaseless PTFE-valve. By using a multi-way valve the vacuum in the collecting tube can be released and the sample removed. After replacing the collecting tube the unit can be again evacuated without interfering with the vacuum in the column.

### Spider Collector

A simple device for the collection of several fractions. Change of the tubes is done by turning the lower part of the collector. Various numbers and sizes of collecting vessels can be supplied. Particular useful when working with substances of low evaporating temperatures as the collecting vessels can be dipped into a Dewar flask. It is also possible to take sample with a syringe from the drop-spit

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



### Low Temperature SPALTROHR<sup>®</sup>-column HMS 1000 T

Separation of mixtures by distillation of low boiling hydrocarbons and other gases, e. g. methane to heptanes, olefins, mixtures of saturated and unsaturated hydrocarbons, can easily be effected with SPALTROHR<sup>®</sup>-columns in laboratory scale in the petrochemical industry and for scientific researches. SPALTROHR<sup>®</sup>-low temperature distillation units enable on one side to run comparing distillations for technical processes in the laboratory and on the other side also gases can be gained in the laboratory. Furthermore, distillations for analytical purposes are possible as well.

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